



Specification

1.06370.0000 Sodium dihydrogen phosphate anhydrous 99.99 Suprapur®

Specification

Purity (metallic)	≥ 99.99	%
Assay (acidimetric)	≥ 99.0	%
Spec. Absorbance A at 260 nm (2 cm; 10 %)	≤ 0.020	
Chloride (Cl)	≤ 5	ppm
Sulfate (SO ₄)	≤ 30	ppm
Total nitrogen (N)	≤ 10	ppm
As (Arsenic)	≤ 0.05	ppm
Ba (Barium)	≤ 10.0	ppm
Ca (Calcium)	≤ 0.50	ppm
Cd (Cadmium)	≤ 0.010	ppm
Ce (Cerium)	≤ 0.010	ppm
Co (Cobalt)	≤ 0.010	ppm
Cs (Cesium)	≤ 5	ppm
Cu (Copper)	≤ 0.010	ppm
Eu (Europium)	≤ 0.010	ppm
Fe (Iron)	≤ 0.05	ppm
Hg (Mercury)	≤ 0.05	ppm
K (Potassium)	≤ 5.0	ppm
La (Lanthanum)	≤ 0.010	ppm
Li (Lithium)	≤ 1.0	ppm
Mg (Magnesium)	≤ 0.50	ppm
Mn (Manganese)	≤ 0.010	ppm
Ni (Nickel)	≤ 0.010	ppm
Pb (Lead)	≤ 0.010	ppm
Rb (Rubidium)	≤ 5.0	ppm
Sm (Samarium)	≤ 0.010	ppm
Sr (Strontium)	≤ 0.50	ppm
Tl (Thallium)	≤ 0.01	ppm
Y (Yttrium)	≤ 0.010	ppm
Yb (Ytterbium)	≤ 0.010	ppm
Zn (Zinc)	≤ 0.010	ppm

Evelyn Allmann

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.