



# Specification

---

1.06549.0000 Sodium nitrite for analysis EMSURE® ACS, Reag. Ph Eur

---

	Spec. Values	
Assay (manganometric)	≥ 99.0	%
Insoluble matter	≤ 0.01	%
Chloride (Cl)	≤ 0.005	%
Sulfate (SO <sub>4</sub> )	≤ 0.005	%
Heavy metals (as Pb)	≤ 0.001	%
Ca (Calcium)	≤ 0.002	%
Fe (Iron)	≤ 0.001	%
K (Potassium)	≤ 0.001	%

Corresponds to ACS, Reag. Ph Eur

Kristin Haufschildt

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.