



Specification

8.56158.0000 Fmoc-D-Glu(OtBu)-Wang resin Novabiochem®

Specification		
Color (visual)	white to yellow to beige	
Appearance of substance (visual)	beads	
Loading (Photometric determination of the Fmoc-chromophore liberated upon treatment with DBU/DMF)	0.40 - 1.00	mmol/g
Swelling Volume (in DMF)	lot specific result	ml/g
Identity (of the substitution)	passes test	

The polymer matrix is copoly (styrene-1% DVB), 200 - 400 mesh. The resin is endcapped after coupling the amino acid derivative.

Dr. Jörg Bauer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.