

# Product Specifications

Product name:

**Alkaline Phosphatase-Anti-Alkaline Phosphatase antibody produced in mouse**

monoclonal AP1B9, purified immunoglobulin, buffered aqueous solution

Product no: **A7827**

## Test:

## Specifications:

### Country of origin

Antibody Origin

Israel

Alkaline Phosphatase Origin

South Africa

BSA Origin

Australia

### Appearance

Form

Liquid

Clarity

Clear

Slightly hazy

Color

Light yellow

Slightly yellowish

Odor

None

### Immunohistochemistry

Tissue

Human appendix

Human tonsil

Fixation

Formalin fixed paraffin embedded

Primary antibody

A2547- Mc anti alpha Smooth Muscle Actin

Bridging antibody

M6898- APA to Mouse Fab ads/w HIgG (Gt)

Substrate

Fast Red

Working dilution

1/40 - 1/80

### Recommended Retest

Recommended Retest period

3 Years