

Product Specifications

Product name:

Anti-Glutathione-S-Transferase (GST)-Agarose antibody produced in rabbit

IgG fraction of antiserum, PBS suspension

Product no: **A8580**

Test:

Specifications:

Country of origin

Antibody origin

Israel

Appearance

Form

Suspension

Clarity

Suspension

Color

White

Odor

None

IgG concentration/ml resin

IgG concentration / 1 ml resin

7 - 8 mg

Binding Capacity

Antigen/Cell extract/Lysate

E. Coli expressing GST tagged protein

Binding capacity

at least 100 ug/ml

Recommended Retest

Recommended Retest period

3 Years