

Product Specifications

Product name:

Anti-Actin, Cardiac antibody, Mouse monoclonal
clone AC1-20.4.2, purified from hybridoma cell culture

Product no: **A9357**

Test:

Appearance

Form
Clarity
Color
Odor

Specifications:

Liquid
Clear
Colorless
None

Protein Concentration by OD at 280nm

IgG concentration 1.0 - 1.5 mg/ml

Immunohistochemistry

tissue
fixation
enzyme
Biotinylated secondary antibody

Rat heart
Formalin fixed paraffin embedded
0.1% pronase
B7151 Anti Mouse IgG (Fab)-Biotin
B0529 Anti Mouse IgG(Fab)-Biotin
E2886 ExtrAvidin Peroxidase
AEC
10 - 20 µg/ml

ExtrAvidin
substrate

working concentration

Recommended Retest

Recommended Retest period

4 Years