

Product Specifications

Product name:

Anti-Calbindin-D-28K (KD-15) antibody produced in rabbit

affinity isolated antibody, buffered aqueous solution

Product no: **C7354**

Test:

Specifications:

Country of origin

Antibody origin

Israel

BSA Origin

Australia

Appearance

Form

Liquid

Clarity

Clear

Color

Colorless

Protein Concentration

Antibody concentration

0.3 - 0.6 mg/ml

Immunoblotting

Antigen

rat recombinant calbindin D-28K(Swant).

Second antibody

A9919 Anti Rabbit IgG(WM)-Alk.Phos.

Substrate

NBT/ BCIP

working dilution

1/1000 - 1/2000

vs human calretinin

Negative

Immunohistochemistry

Tissue

Rat cerebellum

Human brain

Chicken cerebellum

Fixation

Formalin fixed paraffin embedded

Enzyme

0.1% Trypsin

Second antibody

Anti Rabbit IgG (WM)-Biotin

B8895 Anti Rabbit IgG (WM)-Biotin

ExtrAvidin

E2886 ExtrAvidin Peroxidase

Substrate

AEC

Working dilution

1/500 - 1/1000

Recommended Retest

Recommended Retest period

4 Years