

Product Specifications

Product name:

Monoclonal Anti-CD14 antibody produced in mouse

monoclonal UCHM-1, purified immunoglobulin, buffered aqueous solution

Product no: **C7673**

Test:

Specifications:

Country of origin

Antibody origin

Israel

BSA Origin

Australia

Appearance

Form

Liquid

Clarity

Clear

Color

Colorless

Odor

None

Protein Concentration by OD at 280nm

Immunoglobulin content

50 - 200 µg/ml

FACScan Analysis-Indirect

Cell source

Human blood

Second antibody

Anti Mouse IgG (Fab) FACS grade -FITC
F2883 Anti M IgG (WM) ads/Hu Prot.- FITC
F2772 Anti Mouse IgG (Fc) F(ab')₂ -FITC

Sample volume/test

Recommended Retest

Recommended Retest period

3 Years