

Product Specifications

Product name:

Monoclonal Anti-Connexin-43 antibody produced in mouse

monoclonal CXN-6, ascites fluid

Product no: **C8093**

Test:

Specifications:

Country of origin

Antibody origin

Israel

Appearance

Form

Liquid

Clarity

Clear

Slightly hazy

Color

Yellow brown

Odor

None

Protein concentration by Biuret

Protein concentration by Biuret

15 - 55 mg/ml

Isotype concentration by RID

Mouse IgM concentration

0.8 - 1.6 mg/ml

Immunoblotting

Cell tissue extract/Antigen

Mouse brain extract

Second antibody

A9917 Anti Mouse IgG (Fab)-Peroxidase

A2304 Anti Mouse IgG (Fab)-Peroxidase

Substrate

ECL

Working dilution

1/8000 - 1/16000

MW detected (specific band) approx.

43 kDa

Recommended Retest

Recommended Retest period

5 Years