

# Product Specifications

Product name:

**Anti-EDEM1 antibody produced in rabbit**

affinity isolated antibody, buffered aqueous solution

Product no: **E8406**

## Test:

## Specifications:

### Country of origin

Antibody origin Israel

### Appearance

Form Liquid  
Clarity Clear  
Color Colorless  
Odor None

### Protein Concentration by OD at 280nm

IgG concentration 1 - 1.5 mg/ml

### Immunoblotting (chemiluminescent)

Cells 293T cells expressing human EDEM1  
Second antibody A0545 Anti Rabbit IgG-Peroxidase  
Working concentration 0.5 - 1 µg/ml  
MW detected (specific band) approx. ~75kDa  
~75kDa and ~65 kDa

### Immunoprecipitation Assay

Agarose P3476 Protein A agarose  
Extract/lysate 293-T cells expressing mouse EDEM1  
Antibody precipitation amount 0.5 - 1 µg

### Recommended Retest

Recommended Retest period 4 Years