

Product Specifications

Product name:

Anti-Fibroblast Growth Factor Receptor-3, Cytoplasmic antibody produced in rabbit

affinity isolated antibody, buffered aqueous solution

Product no: **F0425**

Test: Specifications:

Country of origin

Antibody origin Israel
BSA Origin Australia

Appearance

Form Liquid
Clarity Clear
Color Colorless
Odor None

Protein Concentration

IgG concentration 1.0 - 1.5 mg/ml

Immunoblotting

Cell/tissue extract 293-T human embryonic kidney cells
Second antibody A9919 Anti Rabbit IgG(WM)-Alk.Phos.
Substrate NBT/ BCIP
Working dilution 1/1000 - 1/2000
MW detected (specific band) 110 - 120 kDa

Immunohistochemistry

Tissue Human umbilical cord
Human placenta
Fixation Formalin fixed paraffin embedded
Enzyme 0.1% Trypsin
Biotinylated second antibody B8895 Anti Rabbit IgG (WM)-Biotin
ExtrAvidin E2886 ExtrAvidin Peroxidase
Substrate AEC
Working dilution 1/250 - 1/500

Recommended Retest

Recommended Retest period 4 Years