

Product Specifications

Product name:

Monoclonal Anti-Green Fluorescent Protein (GFP) antibody produced in mouse

monoclonal GSN149, purified immunoglobulin, buffered aqueous solution

Product no: **G1546**

Test:

Specifications:

Country of origin

Antibody origin

Israel

Appearance

Form

Liquid

Clarity

Clear

Color

Colorless

Odor

None

Protein Concentration by OD at 280nm

IgG concentration

2 - 2.5 mg/ml

Immunoblotting (chemiluminescent)

Cells

293T expressing E-GFP N-terminal

Second antibody

A9917 Anti Mouse IgG (Fab)-Peroxidase

Working concentration

1 - 2 µg/ml

MW detected (specific band) approx.

at~ 65 kDa

Immunoblotting (chemiluminescent)

extract/antigen

Recombinant GFP, purified

Antibody concentration

2 µg/ml

Second antibody code#

A9917 Anti Mouse IgG (Fab)-Peroxidase

Antigen detected

6 - 9 ng/lane

MW detected (specific band)

27 kDa

Recommended Retest

Recommended Retest period

4 Years