

Product Specifications

Product name:

Anti- β -Glucuronidase (C-Terminal) antibody produced in rabbit

~ 1.5mg/mL, affinity isolated antibody, buffered aqueous solution

Product no: **G5545**

Test:

Specifications:

Country of origin

Antibody origin Israel

Appearance

Form Liquid
Clarity Clear
Color Colorless
Odor None

Protein Concentration by OD at 280nm

IgG concentration 1.5 - 2 mg/ml

Immunoblotting (chemiluminescent)

Extract/antigen G7646 Purified beta Glucuronidase(E.Coli)
Second antibody A0545 Anti Rabbit IgG-Peroxidase
Working concentration 1 - 2 μ g/ml
MW detected (specific band) 60 kDa

Recommended Retest

Recommended Retest period 4 Years