

Product Specifications

Product name:

Anti-GRP78/BiP (ET-21) antibody produced in rabbit

IgG fraction of antiserum, buffered aqueous solution

Product no: **G9043**

Test:

Specifications:

Country of origin

Antibody origin

Israel

Appearance

Form

Liquid

Clarity

Clear

Color

Colorless

Odor

None

Protein Concentration by OD at 280nm

IgG concentration

10 - 15 mg/ml

Immunoblotting (chemiluminescent)

Cell extract/antigen

HeLa human cells extract

Second antibody

A0545 Anti Rabbit IgG-Peroxidase

Working dilution

1/3000 - 1/4000

MW detected

band at 78 KDa

band at 78KDa + 73KDa

Immunocytochemistry (IFT)

Cell type

NIH-3T3 mouse fibroblasts

Fixation reagent 1

Methanol

Fixation reagent 2

Acetone

2nd antibody

F9887 Anti Rabbit IgG (WM)-FITC

Working dilution

1/1000 - 1/2000

Recommended Retest

Recommended Retest period

4 Years