

Product Specifications

Product name:

Monoclonal Anti-Neurofilament 160/200 antibody produced in mouse

~ 2mg/mL, monoclonal RMdO20, purified immunoglobulin, buffered aqueous solution

Product no: **N2912**

Test:

Specifications:

Country of origin

Antibody origin

Israel

Immunoglobulin Isotype

Immunoglobulin isotype

Mouse IgG1

Appearance

Form

Liquid

Clarity

Clear

Color

Colorless

Odor

None

Protein Concentration by OD at 280nm

IgG concentration

2 - 2.5 mg/ml

Immunoblotting (chemiluminescent)

extract/antigen

Rat brain extract (S1 fraction)

Rat brain extract

Second antibody code#

A9917 Anti Mouse IgG (Fab)-Peroxidase

A2304 Anti Mouse IgG (Fab)-Peroxidase

ECL-Substrate

ECL

Working concentration

1 - 2 µg/ml

MW detected (specific band)

Bands at 160 KDa and 200 KDa

Recommended Retest

Recommended Retest period

4 Years