

Product Specifications

Product name:

Anti-Sheep IgG (whole molecule) antibody produced in rabbit

whole antiserum

Product no: **S1140**

Test: Specifications:

Country of origin

Antibody origin Israel

Appearance

Form Liquid
Clarity Clear
Color Yellow brown
Odor None

Protein concentration by Biuret

Protein concentration by Biuret 40 - 100 mg/ml

Immunoelectrophoresis (IEP)-Specificity

vs sheep serum one line
one split line
vs sheep IgG one line
one split line

Antibody concentration by precipitin

Antigen Sheep Serum
Antibody concentration 2.8 - 4.2 mg/ml

ELISA Indirect

Antigen IgG from pooled sheep serum
antigen coating conc. 5 µg/ml
Second antibody A0545 Anti Rabbit IgG-Peroxidase
substrate OPD
wave length 450 nm
Working dilution 1/100000 - 1/150000

Recommended Retest

Recommended Retest period 5 Years