

Product Specifications

Product name:

Anti-DIS3 antibody produced in rabbit
~1.0 mg/mL, affinity isolated antibody

Product no: **SAB4200402**

Test:**Specifications:****Country of origin**

Antibody origin

Israel

Appearance

Form

Liquid

Clarity

Clear

Color

Colorless

Odor

None

Protein Concentration by OD at 280nm

IgG concentration

1.0 - 1.5 mg/ml

Immunoblotting (chemiluminescent)

Cell /Tissue /antigen

293-T cells extract

Second antibody

A0545 Anti Rabbit IgG-Peroxidase

Working concentration

2.5 - 5 µg/ml

MW detected (specif.band) approx

109 kDa

Immunocytochemistry (IFT)

Cell type

HeLa cells (human)

Fixation reagent 1

Methanol

Fixation reagent 2

Acetone/ Methanol 1:1

Second antibody

Anti Rabbit IgG-Cy3

Working concentration

0.5 - 1 µg/ml

Recommended Retest

Recommended Retest period

4 Years

Storage conditions

Storage Temperature:

-20°C