

Product Specifications

Product name:

Monoclonal Anti-Human IgM (μ - chain specific)-Alkaline Phosphatase
clone MB-11, Reagents designed and manufactured under cGMP controls suitable for use in an IVD application.

Product no: **SRE0058**

Test:	Specifications:
Appearance	
Form	Liquid
Clarity	Clear
	Slightly hazy
Color	Faint yellow
	Light yellow
	Yellow
	Dark yellow
Odor	None
Conjugate concentration	
IgG content	1.5 - 2.5 mg/ml
ELISA Direct	
Antigen	Human IgM
Antigen coating concentration	5 μ g/ml
Substrate	pNPP
Titer	1/35000 - 1/45000
Dot Blot	
Antigen	Human IgM
Substrate	NBT/ BCIP
Working dilution	1/20000 - 1/40000
Antigen concentration	2.5 - 20 ng/dot
Dot Blot (chemiluminescent)	
Antigen	Human IgM
Working dilution	1/30000 - 1/60000
Antigen concentration	2.5 - 20 ng/dot
Recommended Retest	
Recommended Retest period	2 Years
Country of origin	
Antibody origin	Israel
Alkaline Phosphatase origin	South Africa