

# Product Specifications

Product name:

**Monoclonal Anti-TBP antibody produced in mouse**

~ 2.0mg/mL, monoclonal 58C9, purified immunoglobulin, buffered aqueous solution

Product no: **T1827**

**Test:**

**Specifications:**

**Country of origin**

Antibody origin Israel

**Appearance**

Form Liquid

Clarity Clear

Color Colorless

Odor None

**Protein Concentration by OD at 280nm**

IgG concentration 2.0 - 2.5 mg/ml

**Immunoblotting (chemiluminescent)**

Cell /Tissue /antigen D.Mel cells nuclear extract

Second antibody A2304 Anti Mouse IgG (Fab)-Peroxidase

A9917 Anti Mouse IgG (Fab)-Peroxidase

Working concentration 1 - 2 µg/ml

MW detected (specific band) approx. 42 kDa

**Recommended Retest**

Recommended Retest period 4 Years