

Product Specifications

Product name:

Monoclonal Anti-Troponin T antibody produced in mouse

monoclonal JLT-12, ascites fluid

Product no: **T6277**

Test:

Specifications:

Country of origin

Antibody origin

Israel

Appearance

Form

Liquid

Clarity

Slightly hazy

Clear

Color

Yellow brown

Odor

None

Protein concentration by Biuret

Protein concentration by Biuret

10 - 55 mg/ml

Isotype concentration by RID

Mouse IgG1 concentration

1 - 3 mg/ml

Immunoblotting

Cell tissue extract/Antigen

Rabbit leg skeletal muscle

Second antibody

A9917 Anti Mouse IgG (Fab)-Peroxidase

A2304 Anti Mouse IgG (Fab)-Peroxidase

Substrate

ECL

Working dilution

1/200 - 1/400

MW detected (specific band) approx.

38 kDa

Recommended Retest

Recommended Retest period

5 Years