

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Prostate Specific Antigen from human semen – buffered aqueous solution

Product Number:

P3235

MDL:

MFCD00165544

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
mg protein/ml (UV)	0.2 - 0.6
Purity by SDS Electrophoresis	≥ 99 %
Tested For Infectious Agents	Tested

Specification: PRD.0.ZQ5.10000093541

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.