

Product Specification

Product Name : Anti-Fibronectin antibody, Mouse monoclonal clone IST-4, purified from hybridoma cell culture
Product Number: F0916
MDL Number : MFCD00162291
Storage Temperature : -20°C

Test	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Protein Concentration(OD280nm)	2.0 - 4.0 mg/mL
WB-Cell Line/Tissue Extract	HS-68 Human Foreskin Fibroblasts Cells
WB-Western Blot Titer	1:500-1:1000
ICC/IF Cells	HS-68 Human Foreskin Fibroblasts Cells
ICC/IF Fixation Reagent 1	Methanol
ICC/IF-Fixation Reagent 2	Acetone
ICC/IF-Titer	1:100-1:200
Recommended Retest Period 5 Years	-----

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1223863

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

