

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Haptoglobin from pooled human plasma – ≥95% (SDS-PAGE), lyophilized powder

Product Number:

H3536

CAS Number:

9087-69-8

MDL:

MFCD00131259

Formula Weight:

86,000 g/mol

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
10 mg/mL, H ₂ O	
% Protein (UV)	≥ 50
Purity by SDS Electrophoresis	≥ 95 %
Tested For Infectious Agents	Tested

Specification: PRD.0.ZQ5.10000104445

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

