

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Protein A from Staphylococcus aureus – Soluble, Cowan strain, recombinant, expressed in E. coli, aqueous solution, ≥95% (HPLC)

Product Number: **P7837**

MDL: MFCD00081934

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Protein by UV (based on E1% at 275 nm = 1.65)	≥ 50 mg/ml
Purity (HPLC)	≥ 95 %
SDS Electrophoresis One major band (44,000 to 46,000 Daltons)	Conforms
Endotoxin Level	≤ 0.5 EU/mg

Specification: PRD.1.ZQ5.10000043720

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.