

Product Specification

Product Name :

Albumin human Exbumin™, Recombinant, lyophilized powder, Animal Component Free, suitable for cell culture, low endotoxin

Product Number: SAE0160**Storage Temperature :** -20°C

Test	Specification
Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder
Solubility (Color)	Very Faint Yellow to Yellow to Brown-Yellow
Solubility (Turbidity) 50 mg/mL, H ₂ O	Clear to Hazy
Water (by Karl Fischer)	≤ 12%
Cell Culture Test	Pass
Endotoxin Level	≤ 1EU/mg
Total Albumin Content by Native-PAGE	≥ 96%

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1191663

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

