

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Hypo-Opticlear Human Sera

Product Number: SAE0012

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Brown and Colorless to Dark Brown-Yellow
Appearance (Form)	Liquid
Cholesterol	≤ 10 mg/dl
Triglyceride	≤ 10 mg/dl
Calcium (Ca) (mg/dl)	7.5 - 10.5
Total Protein g/dl	≥ 6.0
pH	7.0 - 8.5
Chloride (mEq/L)	95 - 110
Potassium (mEq/L)	3.5 - 5.0
Tested For Infectious Agents	Tested
Source	Conforms
Plasma from donors obtained and processed in the United States	
Sodium (Na) (mEq/L)	130 - 165

Specification: PRD.1.ZQ5.10000073526

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.