

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Insulin-like Growth Factor-I human - IGF-I, recombinant, expressed in E. coli, lyophilized powder, suitable for cell culture

Product Number:

I3769

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Form)	Powder
Lyophilized	
Purity (HPLC)	≥ 95.0 %
Sum the 2 peaks in case the material is oxygenated.	
Activity	0.100 - 10.000 ng/ml
ED50	
Endotoxin Level	≤ 1.00 EU/ug
Recommended Retest Period	-----
3 Years	

Specification: PRE.1.ZQ5.10000021190

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

