

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Hydroxyphenyl fluorescein solution - 5 mM in DMF

Product Number: **H4290**

Formula: C<sub>26</sub>H<sub>16</sub>O<sub>6</sub>  
Formula Weight: 424.40 g/mol  
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Yellow to Very Dark Yellow and Light Orange to Orange
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
Fluorescence	512 - 518 nm
Emission Maximum	
Fluorescence	480 - 493 nm
Excitation Maximum	
Purity (HPLC)	≥ 98 %

Specification: PRD.2.ZQ5.10000036219

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.