

Product Specification

Product Name : Serum minus IgG human lyophilized powder
Product Number: S5143
Storage Temperature : +2°C to +8°C

Test	Specification
Appearance (Form)	Powder
Appearance (Color)	Off White
Solubility (Solvent)	Water
Solubility (Volume)	0.5 - 0.5 ml
Solubility (Color)	Very Light Brown to Very Dark Brown and Dark Brown-Red to Very Dark Brown-Red
Solubility (Turbidity)	Clear to Hazy
Protein by Biuret	50 - 100 mg/mL
Human Source Testing (Source)	-----
Human Source (HIV)	Negative
Human Source (HbSAg)	Negative
Human Source (HCV)	Negative
Human Source (Syphilis)	Negative
Recommended Retest Period 5 Years	-----

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1223321

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

