

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Biuret reagent – 150–1,000 µg/mL protein detection level

Product Number:

**B3934**

MDL:

MFCD00007946

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Turbidity)

Clear

Appearance (Color)

Conforms

Pale Blue

Appearance (Form)

Liquid

Suitability

Suitable

Tested and found suitable for use in the quantitation of protein.

Specification: PRD.1.ZQ5.10000006164

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.