

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:


L- α -Phosphatidylcholine - from egg yolk, Type XVI-E, $\geq 99\%$ (TLC), lyophilized powder

Product Number:

P3556

CAS Number:

8002-43-5



Storage Temperature:

-20 °C

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
100 mg/mL, CHCl ₃	
Phosphorus (P)	3.4 - 4.4 %
Purity (TLC)	$\geq 99\%$
Note	-----

Natural product, the composition of palmitic, oleic, linoleic and stearic acid will vary.

Specification: PRD.0.ZQ5.10000040060

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.