

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
IgG from human serum – technical grade, ≥80% (SDS–PAGE), buffered aqueous solution

Product Number: I8640
MDL: MFCD00163923
Storage Temperature: -20 °C

TEST	Specification
mg protein/ml (UV)	10.0 - 14.0
ID by Immuno-Electrophoresis	Pass
Single major arc of precipitation in the gamma region versus anti-human whole serum and anti-human IgG. Minor gamma arcs may be present versus anti-human whole serum	
Purity by SDS-PAGE	≥ 80 %
Tested For Infectious Agents	Tested

Specification: PRD.0.ZQ5.10000049413

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.