

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Trimethylolpropane trimethacrylate – contains 250 ppm monomethyl ether hydroquinone as inhibitor, technical grade

Product Number:

246840

CAS Number:

3290-92-4

MDL:


MFCD00008588

Formula:

C₁₈H₂₆O₆

Formula Weight:

338.40 g/mol



TEST

Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Viscous Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

Typically 90%

MEHQ as Inhibitor

80 - 350 ppm

Specification: PRD.1.ZQ5.10000040145

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

