

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Leukocyte Alkaline Phosphatase Kit – based on naphthol AS-BI and fast blue BB salt

Product Number:

86C

Storage Temperature:

2 - 8 °C

TEST

Specification

Certificate of Analysis

Confirmed

The individual components of this kit conforms to the MilliporeSigma quality testing and performance specifications as outlined in the package insert

Specification: PRD.0.ZQ5.10000103736

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

