

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Reversible Protein Detection Kit for Membranes and Polyacrylamide Gels

Product Number:

RPROB

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance

Purple when 10 ml Solution A is mixed with 10 ml Solution B

Conforms

Specification: PRD.0.ZQ5.10000093687

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.