

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
OptiPrep™ Density Gradient Medium - used for cell and subcellular organelle isolation

Product Number: **D1556**

MDL: MFCD00867965



TEST	Specification
Bulk Density	1.319 - 1.321 g/ml
Endotoxin Level	≤ 1 EU/ml
Sterility	Pass
Trademark	Confirmed

OptiPrep is a trademark of Serumwerk Bernburg AG

Specification: PRD.2.ZQ5.10000033831

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

