

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Interferon  $\beta$  from mouse –  $\geq 95\%$ , recombinant, expressed in E. coli, buffered aqueous solution, cell culture testedProduct Number: **I9032**

MDL: MFCD02189136

Storage Temperature: -70 °C

TEST	Specification
Vial Content	Conforms
Not less than 100,000 units/vial	
Specific Activity	
Report Result (units/mg)	
Appearance	Conforms
Solution in 20 mM HEPES, pH 6.0, 0.5 M NaCl, 6% glycerol, and 0.1% bovine serum albumin.	
Purity	Conforms
Not less than 95%	
Concentration	
Report Result (units/ml)	

Specification: PRD.0.ZQ5.10000074294

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.