


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Piperazine diacrylamide - For acrylamide gel electrophoresis

Product Number: **D1538**
CAS Number: 6342-17-2
MDL: MFCD00077661
Formula: C₁₀H₁₄N₂O₂
Formula Weight: 194.23 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear to Hazy
at 40 mg plus 4 mL of H ₂ O	
Suitability	Suitable
Suitable for use in polyacrylamide gel electrophoresis	
Purity (HPLC)	≥ 99 %

Specification: PRD.1.ZQ5.10000040991

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.