

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:


Tris Acetate-EDTA buffer – 10× concentrate, BioReagent, for molecular biology, non-sterile; 0.2 µm filtered

Product Number:

T9650

MDL:

MFCD00236357



TEST	Specification
------	---------------

Appearance (Turbidity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Solution

pH

8.2 - 8.4

at 25 degrees Celcius

RNAse Detection

None Detected

Protease Detection

None Detected

Suitability

Suitable

for use in Gel Electrophoresis after dilution to working concentration, and no DNase/Nickase activity detected

Recommended Retest Period

Complete

1 year

Specification: PRD.2.ZQ5.10000070126

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.