

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Rabbit IgG (whole molecule)-Biotin antibody produced in goat - affinity isolated antibody, buffered aqueous solution

Product Number:

B7389

MDL:

MFCD00162783

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Turbidity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Solution

Protein Content

0.1 - 2.0 mg/ml

ELISA Titer

≥ 20000

Minimum 1:20,000 working dilution versus 5 Åµg/ml normal rabbit IgG.

Specification: PRD.1.ZQ5.10000068951

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

